

CERAMAT

Technical Data Sheet

Product Name: Phase Pure Hydroxyapatite

Mol. Formula: $[\text{Ca}_{10}(\text{PO}_4)_6(\text{OH})_2]$

Mol. Weight: 1004.62 gm/mol

CAS No.: 12167-74-7

Sl. No.	Tests	Tests methods	Specifications
1.	Phase Analysis	XRD	Conforms with ICDD file no. 09-0432
2.	Appearance	Visually	Milky white with no other-coloured impurity.
3.	Ca/P ratio	ICP-MS/AES/OES	$1.663 \leq \text{Ca/P} \leq 1.714$
4.	Trace elements	ICP-MS/AES/OES	As \leq 3 ppm Cd \leq 5 ppm Hg \leq 5 ppm Pb \leq 30 ppm
5.	Particle Size	Sieve Analysis	$<75\mu\text{m}$
6.	Crystallite Size	XRD Debye- Scherrer	15-45 nm
7.	Heavy metals	ICP-MS/AES/OES	Total heavy metals \leq 30 ppm
8.	Other elements	ICP-MS/AES/OES	Fe \leq 100 ppm
9.	Bulk Density	Manually (Mass/Volume)	0.2 to 0.8 gm/cc

Ceramat Private Limited (CIN: U26990MH2021PTC370837, GSTN: 27AAJCC8882E1Z9)
Office Address: Unit No. 14, Sethia Industrial Park, S.No. 39/1, Bilalpada, Vasai East-401208,
Maharashtra, India.

Tel: +91 8857861496; Web site: www.tslceramat.com ;Email- info@tsl-ceramat.com

CMPL/L3/C.Comm/06 Ver 01 Date:19-06-2023